




**THUASNE**

# Epi PROMASTER & Epi-Med PROMASTER



### Anatomical fit

- Supports the joint area
- Prevents slipping

### Soft knit fabric

- Maximum comfort in the elbow crease



### FlexAir pads

- Relieve pressure on the tendon insertions in the forearm muscles



## Epi PROMASTER

Product Number **35 080**

Colour	Version	Dimensions		Qty.
		Size	Forearm circumference	
Dark grey/Blue	Universal	XS	21-23 cm	
		S	23-25 cm	
		M	25-27 cm	
		L	27-29 cm	
		XL	29-31 cm	
		XXL	31-34 cm	

## Epi-Med PROMASTER

Product Number **35 087**

Colour	Version	Dimensions		Qty.
		Size	Forearm circumference	
Dark grey/Blue	Universal	XS	22-23 cm	
		S	24-25 cm	
		M	26-27 cm	
		L	28-30 cm	
		XL	31-34 cm	

**BMI**  
 BERRETTA MEDICAL INC.

The exclusive Canadian distributor for Thuasne USA

545 St. Jacques O. Montreal, Qc H4A 2E3  
 O: 514.369.3311

Toll Free / Sans Frais : 1.877.222.3311 ✶ Fax Toll Free / Telec Sans Frais : 1.877.527.1911  
 www.bmiortho.com