



THUASNE

Genu PROMASTER



Y pad

- Patella support
- Surrounding tissue is massaged, extending right up to the thighs
- The distal area of the pad is cut out to relieve pressure on the edge of the shinbone



Comfortable, firm fit

- Gradual change in compression prevents pinching
- Breathable knit fabric with 3-D texture for comfort

BMI
BERRETTA MEDICAL INC.

The exclusive Canadian distributor for Thuasne USA



Grips on both sides

- Make it easier to apply the soft brace
- Patent pending



Welded spiral springs

- Enhanced stabilisation of the knee joint
- Gently support knee movement, without obstructing it
- Stop the brace from slipping

Genu PROMASTER
Product Number **35 040**

Colour	Version	Dimensions			Qty.
		Size	Thigh circumference	Calf circumference	
Dark grey/Blue	Universal	XXS	35-38 cm	25-28 cm	
		XS	38-41 cm	28-31 cm	
		S	41-44 cm	31-34 cm	
		M	44-47 cm	34-37 cm	
		M+	47-50 cm	34-37 cm	
		L	47-50 cm	37-40 cm	
		L+	50-53 cm	37-40 cm	
		XL	50-53 cm	40-43 cm	
		XL+	53-56 cm	40-43 cm	
		XXL	53-56 cm	43-46 cm	
XXL+	56-59 cm	43-46 cm			

5545 St. Jacques O. Montreal, Qc H4A 2E3
O: 514.369.3311

Toll Free / Sans Frais : 1.877.222.3311 ☘ Fax Toll Free / Telec Sans Frais : 1.877.527.1911
www.bmiortho.com

