





THUASNE

Malleo PROMASTER & Malleo PROMASTER Aktiv



Anatomically shaped pads

- Effectively support the ankle
- FlexAir structure guarantees the circulation of air

Patented pronation strap Malleo PROMASTER Aktiv

- Enhanced forefoot stability
- Additional protection against the foot buckling inwards
- Thus preventing resulting additional stress on the ankle ligaments and tendons

Malleo PROMASTER Product Number 35 020

Colour	Version	Dimensions		Qty. RT	Qty. LT
		Size	Ankle circumference		
Dark grey/Blue	Right or left	XS	17-19 cm		
		S	19-21 cm		
		M	21-23 cm		
		L	23-25 cm		
		XL	25-27 cm		
		XXL	27-29 cm		

Malleo PROMASTER Aktiv Product Number 35 022

Colour	Version	Dimensions		Qty. RT	Qty. LT
		Size	Ankle circumference		
Dark grey/Blue	Right or left	XS	17-19 cm		
		S	19-21 cm		
		M	21-23 cm		
		L	23-25 cm		
		XL	25-27 cm		
		XXL	27-29 cm		



BMI
BERRETTA MEDICAL INC.

The exclusive Canadian distributor for Thuasne USA

545 St. Jacques O. Montreal, Qc H4A 2E3
O: 514.369.3311

Toll Free / Sans Frais : 1.877.222.3311 Fax Toll Free / Telec Sans Frais : 1.877.527.1911
www.bmiortho.com